Lifecosm COVID-19 Antigen Test Cassette

Instruction for use
Gebrauchsanweisung
Mode d'emploi
Instrucciones de uso
Istruzioni per uso
Gebruiksaanwijzing

REF:RC-HM01
Version No.:4.1
Effective Date:
August 1, 2022

Lifecosm Biotech Limited

Add: 3B, Fusen Building, 669 Shuangming Avenue, Guangming District, Shenzhen, China.

Email: info@lifecosm.com

SUNGO Europe B.V.

Add: Olympisch Stadion 24, 1076DE Amsterdam, Netherlands

Contact: Yana Zhang (Ms.)

Tel/Fax: +31 (0) 2021 11106

E-mail: yan.zhang@sungoglobal.com

ENGLISH

[INTENDED USE]

The COVID-19 Antigen Test Cassette is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasopharyngeal swab, oropharyngeal swab, nasal swab, or saliva from individuals who are suspected of COVID-19 by their healthcare provider. Results are for the identification of SARS-CoV-2 nucleocapsid antigen. Antigen is generally detectable in oropharyngeal swab, nasal swab, or saliva during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease. Negative results do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results should be considered in the context of a patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19, and confirmed with a molecular assay, if necessary for patient management. The COVID-19 Antigen Rapid Test Cassette intended for use by medical professionals or trained operators who are proficient in performing lateral flow tests. The product may be used in any laboratory and non-laboratory environment that meets the requirements specified in the Instructions for Use and local regulation.

[WARNINGS AND PRECAUTIONS]

- For in vitro diagnostic use only.
- For healthcare professionals and individuals trained in point of care settings.
- Do not use this product as the sole basis to diagnose or exclude SARS-CoV-2 infection or to inform infection status of COVID-19.
- Do not use this product after the expiration date.
- Please read all the information in this leaflet before performing the test.
- The test cassette should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test cassette should be discarded according to federal, state and local regulations.

[COMPOSITION]

Materials Provided

- 25 Test Cassettes: each cassette with desiccant in individual foil pouch
- 25 Sterilized Swabs: single use swab for specimen collection
- 25 Extraction Tubes: containing 0.4mL of extraction reagent
- 25 Dropper Tips
- 1 Work Station
- 1 Package Insert

Materials Required but not Provided: Timer


[STORAGE AND STABILITY]

- Store as packaged in the sealed pouch at temperature (4-30°C or 40-86°F). The kit is stable within the expiration date printed on the labeling.
- Once open the pouch, the test should be used within one hour. Prolonged exposure to hot and humid environment will cause product deterioration.
- The LOT and the expiration date were printed on the labeling.

[SPECIMEN]

Specimens obtained during early symptom onset contain the highest viral titers; specimens obtained after five days of symptoms are more likely to produce negative results when compared to an RT-PCR assay. Inadequate specimen collection, improper specimen handling and/or transport may yield false results; therefore, training in specimen collection is highly recommended due to the importance of specimen quality to obtain accurate test results. Acceptable specimen type for testing is a direct swab specimen or a swab in viral transport media (VTM) without denaturing agents. Use freshly collected direct swab specimens for best test performance. Prepare the extraction tube according to the Test Procedure and use the sterile swab provided in the kit for specimen collection.

Nasopharyngeal Swab Specimen Collection


- Remove the swab from the package.
- Tilt patient's head back about 70°.
- Insert the swab through the nostril parallel to the palate (not upwards) until resistance is encountered or the distance is equivalent to that from the ear to the nostril of the patient, indicating contact with the nasopharynx. (Swab should reach depth equal to distance from nostrils to outer opening of the ear.) Gently rub and roll the swab. Leave swab in place for several seconds to absorb secretions.
- Slowly remove swab while rotating it. Specimens can be collected from both sides using the same swab, but it is not necessary to collect specimens from both sides if the tip of swab is saturated with fluid from the first collection. If a deviated septum or blockage creates difficulty in obtaining the specimen from one nostril, use the same swab to obtain the specimen from the other nostril.

Oropharyngeal Swab Specimen Collection

Insert swab into the posterior pharynx and tonsillar areas. Rub swab over both tonsillar pillars and posterior oropharynx and avoid touching the tongue, teeth, and gums.

Nasal Swab Specimen Collection

- Take out the swab from the package.
- Tilt patient's head back about 70°.
- While gently rotating the swab, insert swab about 2.5 cm (1 inch) into nostril until resistance is met at turbinates.
- Rotate the swab several times against nasal wall and repeat in another nostril using the same swab.

Saliva Specimen Collection


- Make sure no eating, drinking, smoking, gum chewing 30 minutes before saliva collection.
- Take out the swab from the package.
- Place the swab head between upper and lower molar teeth for 10 seconds to absorb saliva.

Specimen Transport and Storage

Do not return the swab to the original swab packaging. Freshly collected specimens should be processed as soon as possible, but no later than one hour after specimen collection. Specimen collected may be stored at 2-8°C for no more than 24 hours; Store at -70°C for a long time, but avoid repeated freeze-thaw cycles.

[TEST PROCEDURE]

- Note: Allow the test cassettes, reagents and specimens to equilibrate to room temperature (15-30°C or 59-86°F) prior to testing.
- Peel off aluminum foil seal from the top of the extraction tube containing the extraction tube.
 - Put an extraction tube on the work station.
 - Sampling refers to section 'Specimen Collection'.
 - Insert the swab specimen into the extraction tube which contains extraction reagent. Roll the swab at least 5 times while pressing the head against the bottom and side of the extraction tube. Leave the swab in the extraction tube for one minute.
 - Remove the swab while squeezing the sides of the tube to extract the liquid from the swab. The extracted solution will be used as test sample.
 - Cover the extraction tube with a dropper tip tightly.
 - Remove the test cassette from the sealed pouch.
 - Reverse the specimen extraction tube, holding the tube upright, transfer 3 drops (approximately 100 µL) slowly to the specimen well (S) of the test cassette, then start the timer.
 - Wait for colored lines to appear. Interpret the test results at 15 minutes. Do not read results after 20 minutes.

[INTERPRETATION OF RESULTS]

Positive		Two lines appear. One colored line appears at the control region (C), and another colored line appears at the test region (T), regardless of the intensity of the test line.
Negative		One colored line appears at the control region (C), and no line appears at the test region (T).
Invalid		Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test cassette. If the problem persists, discontinue using the lot immediately and contact your local distributor.

[QUALITY CONTROL]

A procedural control is included in the test. A colored line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

[LIMITATIONS]

- The product is limited to provide a qualitative detection. The intensity of the test line does not necessarily correlate to the concentration of the antigen of the specimens.
- Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for patient management decisions.
- A physician must interpret the results in conjunction with the patient's history, physical findings, and other diagnostic procedures.
- A negative result can occur if the quantity of SARS-CoV-2 antigens present in the specimen is below the detection threshold of the assay, or the virus has undergone minor amino acid mutation(s) in the target epitope region recognized by the monoclonal antibodies utilized in the test.

[PERFORMANCE CHARACTERISTICS]

Clinical Performance

For nasopharyngeal swab:

The clinical performance of COVID-19 Antigen Rapid Test Cassette was established in prospective studies with nasopharyngeal swabs collected from 405 individual symptomatic patients (within 7 days of onset) and asymptomatic patients who were suspected of COVID-19. Summary data of COVID-19 Antigen Rapid Test as below:

COVID-19 Antigen	RT-PCR		Total
	Positive	Negative	
LIFECOSM	Positive	158	159
	Negative	6	319
Total	164	314	478

PPA: 96.34% (158/164), (95%CI: 92.2%~98.7%)
NPA: 99.68% (313/314), (95%CI: 98.2%~99.9%)

For nasal swab:

The clinical performance of COVID-19 Antigen Rapid Test Cassette was established in prospective studies with nasal swabs collected from 605 individual symptomatic patients (within 7 days of onset) and asymptomatic patients who were suspected of COVID-19. Summary data of COVID-19 Antigen Rapid Test as below:

COVID-19 Antigen	RT-PCR		Total
	Positive	Negative	
LIFECOSM	Positive	131	134
	Negative	9	462
Total	140	465	605

PPA: 93.6% (131/140), (95%CI: 88.2%~96.6%)
NPA: 99.4% (462/465), (95%CI: 98.1%~99.8%)

For saliva swab:

The clinical performance of COVID-19 Antigen Rapid Test Cassette was established in prospective studies with saliva swabs collected from 585 individual symptomatic patients (within 7 days of onset) and asymptomatic patients who were suspected of COVID-19. Summary data of COVID-19 Antigen Rapid Test as below:

COVID-19 Antigen	RT-PCR		Total
	Positive	Negative	
LIFECOSM	Positive	91	93
	Negative	9	483
Total	100	485	585

PPA: 91% (91/100), (95%CI: 83.8%~95.2%)
NPA: 99.6% (483/485), (95%CI: 98.5%~99.9%)

Limit of Detection (Analytical Sensitivity)

Minimal detection limits of COVID-19 Antigen Test Cassette is 50pg/ml for recombinant COVID-19 protein.

Cross Reactivity (Analytical Specificity)

Cross reactivity was evaluated by testing 32 commensal and pathogenic microorganisms that may be present in the nasal cavity. No cross-reactivity was observed with recombinant MERS-CoV NP protein when tested at the concentration of 50 µg/mL. No cross-reactivity was observed with the following viruses when tested at the concentration of 1.0x10⁶ PFU/mL: Influenza A (H1N1), Influenza A (H1N1pdm09), Influenza A (H3N2), Influenza B (Yamagata), Influenza B (Victoria), Adenovirus (type 1, 2, 3, 5, 7, 55), Human metapneumovirus, Parainfluenza virus (type 1, 2, 3, 4), Respiratory syncytial virus, Enterovirus, Rhinovirus, Human coronavirus 229E, Human coronavirus OC43, Human coronavirus NL63, Human coronavirus HKU1. No cross-reactivity was observed with the following bacteria when tested at the concentration of 1.0x10⁷ CFU/mL: Mycoplasma pneumoniae, Chlamydia pneumoniae, Legionella pneumophila, Haemophilus influenzae, Streptococcus pyogenes (group A), Streptococcus pneumoniae, Candida albicans, Staphylococcus aureus.

Interference

The following potential interference substances were evaluated with the COVID-19 Antigen Rapid Test Cassette at the concentrations listed below and were found not to affect test performance.

Substance	Concentration	Substance	Concentration	Substance	Concentration
Mucin	2%	Arbidol	5 mg/mL	Mupirocin	10 mg/mL
Benzoicaine	5 mg/mL	Phenylephrine	15%	Zanamivir	5 mg/mL
Saline nasal spray	15%	Triamcinolone	10 mg/mL	Ribavirin	5 mg/mL
Oxymetazoline	15%	Whole blood	4%	Dexamethasone	5 mg/mL
Tobramycin	5 µg/mL	Menthol	10 mg/mL	Histamine	10 mg/mL
Oseltamivir phosphate	10 mg/mL	Fluticasone propionate	5%	dihydrochloride	

High-dose Hook Effect

The COVID-19 Antigen Rapid Test Cassette was tested by recombinant COVID-19 protein solutions up to 100µg/ml, and no high-dose hook effect was observed.

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